

ilShin[®]

Ice Flake Machine



The only ice machine specifically
designed for the laboratory

ilShin[®] Europe

iShin[®] Icemaker

The iShin Icemaker is the only Icemaker which is specifically designed for the laboratory. Compared with regular commercial machines the iShin offers some unique design specifications. The cabinet is made of indestructible HDPE (High Density Poly Ethylene). This material will withstand decennia of continues operation!

Features and Benefits



Safety System for ice production

A robust temperature sensor will automatically switch off the ice production as soon as the ice reaches the top of the bin.



WhisperFreeze Technology[®]

Powerful compressor achieves fast ice production with the lowest noise level.



Built-in Water Filter

Make sure your samples are covered in clean and pure flake ice.

Specifications IF300-150	
External Dimensions	700 x 700 x 975 mm (W x D x H)
Bin Storage Dimensions	600 x 300 x 500 mm (W x D x H)
Ice Production	120 kg / 24 hour
Storage capacity	70 kg ice
Refrigeration System	1/3 HP compressor
Refrigerant	R402B, 360g
Electrical Requirements	220-240 V, 50-60 Hz, single phase
Power consumption	0.5kW, 2.3A
Gross Weight	114 kg
Defrost	User friendly drain for defrost water

Your local iShin dealer: