



## Storage Racks with side opening (Uprights)

**Model:**

**DF8510**


**DF8514**

**DF8517**

**DF8520**

**DF8524**


Amount of racks:  
(view from afront)




12 racks




16 racks



16 racks



20 racks



20 racks

Type of racks for 2 inch boxes:



3 deep, 5 high  
Part: UR0066



3 deep, 5 high  
Part: UR0066



4 deep, 5 high  
Part : UR0068



4 deep, 5 high  
Part : UR0068



5 deep, 5 high  
Part : UR0070

Type of racks for 3 inch boxes:



3 deep, 3 high  
Part: UR0067



3 deep, 3 high  
Part: UR0067



4 deep, 3 high  
Part : UR0069



4 deep, 3 high  
Part : UR0069



5 deep, 3 high  
Part : UR0071

Max.of 2 inch:

180 boxes

240 boxes

320 boxes

400 boxes

500 boxes

Max. of 3 inch:


108 boxes

144 boxes

192 boxes

240 boxes

300 boxes



External Rack Dimensions:  
(mm D x W x H)

425x145x283

425x145x283

558x145x283

558x145x283

692x145x283

Internal Rack Dimensions  
(mm D x W x H)

For 2 inch boxes, each opening: 145 x 136 x 51


For 3 inch boxes, each opening: 145 x 136 x 88




## Storage Racks with drawers (Uprights)

**Model:**                      **DF8510**                      **DF8514**                      **DF8517**                      **DF8520**                      **DF8524**


**Amount of racks:**  
(view from afront)




12 racks




16 racks



16 racks



20 racks



20 racks

**Type of racks for 2 inch boxes:**



3 deep, 5 high  
Part: UR1066



3 deep, 5 high  
Part: UR1066



4 deep, 5 high  
Part : UR1068



4 deep, 5 high  
Part : UR1068



5 deep, 5 high  
Part : UR1070

**Type of racks for 3 inch boxes:**



3 deep, 3 high  
Part: UR1067



3 deep, 3 high  
Part: UR1067



4 deep, 3 high  
Part : UR1069



4 deep, 3 high  
Part : UR1069



5 deep, 3 high  
Part : UR1071

**Max.of 2 inch:**

180 boxes

240 boxes

320 boxes

400 boxes

500 boxes

**Max. of 3 inch:**


108 boxes

144 boxes

192 boxes

240 boxes

300 boxes



**External Rack Dimensions:**  
(mm D x W x H)

430x140x290	430x140x290	560x140x290	560x140x290	700x140x290
-------------	-------------	-------------	-------------	-------------

**Internal Rack Dimensions**  
(mm D x W x H)

For 2 inch boxes, each opening: 136 x 136 x 52

For 3 inch boxes, each opening: 136 x 136 x 89



## Storage Racks for Chest Freezers

**Model:**

**DF9007**

**DF9010**

**DF9014**

**DF9017**

**DF9020**

**Amount of racks:**  
(view from above)



12 racks



18 racks



24 racks



30 racks



36 racks

**Type of racks for 2 inch boxes:**



12 high  
Part: CR0062



12 high  
Part: CR0062



12 high  
Part: CR0062



12 high  
Part: CR0062



12 high  
Part: CR0062

**Type of racks for 3 inch boxes:**



8 high  
Part: CR0063



8 high  
Part: CR0063



8 high  
Part: CR0063



8 high  
Part: CR0063



8 high  
Part: CR0063

**Max. of 2 inch:**

144 boxes

216 boxes

288 boxes

360 boxes

432 boxes

**Max. of 3 inch:**


96 boxes

144 boxes

192 boxes

240 boxes

288 boxes



**External Rack Dimensions:**  
(mm D x W x H)


146x146x660	146x146x660	146x146x660	146x146x660	146x146x660
-------------	-------------	-------------	-------------	-------------

**Internal Rack Dimensions:**  
(mm D x W x H)

For 2 inch boxes, each opening: 145 x 138 x 51

For 3 inch boxes, each opening: 137 x 134 x 78

## Storage Boxes for ilShin ULT Freezers



Easy to remove cover

Thick Plastified Cardbox Material

### External dimensions:

External dimensions of 2 inch boxes (model PP0060): 133x133x48 mm (DxWxH)

External dimensions of 3 inch boxes (model PP0061): 133x133x77 mm (DxWxH)

### Internal Dimensions:

Internal dimensions of 2 inch boxes (model PP0060): 124x124x45 mm (DxWxH)

Internal dimensions of 3 inch boxes (model PP0060): 124x124x70 mm (DxWxH)

## Possible dividers:



Dividers are also made of plastified cardbox. Height of all dividers is 35 mm

<b>Model:</b>	<b>Description:</b>
---------------	---------------------

PP1049	Storage box divider 7 x 7, (for vials < 16 mm O.D)
PP1064	Storage box divider 8 x 8, (for vials < 14 mm O.D)
PP1081	Storage box divider 9 x 9, (for vials < 13 mm O.D)
PP1100	Storage box divider 10 x 10, (for vials < 12 mm O.D)